

CINCO VRA PRO 2 System

- 4 Side reinforcers (2 left and 2 right)
- 1 Animated Centering Toy (Same size as side units)
- 1 wireless remote
- 3 power supplies with extensions for each unit
- 2 power supply jumpers for secondary units
- Installation kit
- ASP & ASP 2 Compatible (with primary units)

This is our Flagship VRA system. It allows for control of 5 units with a single wireless remote.



Our VRA PRO System was designed and developed BY an audiologist, FOR audiologists.



Conditioned Play Innovations VRA Systems 2372 W 1200 N Provo, Utah 84601

> Phone: (801) 607-1800 E-mail: sales@vrapro.com Website: www.vrapro.com

Would you like to know more about our VRA Pro System? We can put you in contact with your local sales representative for more information.



Visual Reinforcement Audiometry



VRA PRO 2 SYSTEM OPTIONS

for

testing infants from ages 6 months to 36 months

-Providing Innovative Pediatric Audiology Solutions-



DUAL VRA PRO 2 System

- 2 Side reinforcers (1 left and 1 right)
- 1 Centering Toy/Midline distractor (non-animated)
- 1 wireless remote
- 2 power supplies with extensions for side units
- 25' RCA cable for center toy power
- Installation kit
- ASP & ASP 2 Compatible



TRIO VRA PRO 2 System

- 2 Side reinforcers (1 left and 1 right)
- 1 Animated Centering Toy (Same size as side units)
- 1 wireless remote
- 3 power supplies with extensions for each unit
- Installation kit
- ASP & ASP 2 Compatible



QUAD VRA PRO 2 System

- 4 Side reinforcers (2 left and 2 right)
- 1 Centering Toy/Midline distractor (non-animated)
- 1 wireless remote
- 2 power supplies with extensions for side units
- 2 power supply jumpers for secondary units
- 25' RCA cable for center toy power
- Installation kit
- ASP & ASP 2 Compatible (with primary units)

Features & Advantages of the VRA PRO 2

- Encrypted Bluetooth Wireless Remote No frequency overlap like low frequency RF or need for install of additional sensors like IR systems
- •No Need for Additional Sensors Easier Install
- •No Frequency Overlap-Great for Adjacent Systems/Rooms
- Spacing Saving Design-Sits nearly 3" deeper than similarly sized 3D toy systems
- Expandable system* (Same model only)-Can operate up to 5 units with the same remote

- •Ultra Bright LEDs-Helps keep the child's attention during testing
- Easy Pairing with Remote-Plug it in and push the remote buttons within two minutes
- Centering Box/Mid-Line Distractor Toy Included Easier testing with one audiologist
- Easy Installation-Designed with the installer in mind
- No Interconnecting Cables (between left and right side reinforcers)-Easier Install
- •No Need for Additional Install Kits or Brackets-Standard Pricing